## How to take a Rapid Antigen Test

Most Rapid Antigen Test kits generally follow the instructions below, but not all. Please follow the manufacturer's instructions.





- After 15 minutes look at the square area on the testing tray
- One red line under the C = negative for virus
- Two lines red line under the C and blue line under T = the virus has been detected (even if the blue line is faint)
- One blue line under C = Invalid result
- Report your results on Mycovidrecord.health.nz, or call 0800 222 478 and press option 3
- If your test has detected the virus:
  - Please advise your employer and any other people you may have been in contact with
  - You will need to go into isolation. For up-to-date isolation requirements go to the Advice for people with COVID-19 on the Ministry of Health website.





Not detected















Failed result

Pacific





